

# Performance Evaluation of pH Measurements Made Using Origin™ for Integration with Chest Tube Drainage Systems




Nour Helwa, Diana Khater, Manaswi Sharma

## Introduction


Chest tubes are routinely used following cardiothoracic surgeries<sup>1,2</sup>. **To date, there is no standardized method for chest tube management, particularly when considering criteria for their removal**<sup>1,2</sup>. Prolonged chest tube use increases the risk of iatrogenic infection, length of stay (LOS), pain, and reduces mobility<sup>1,3</sup>.

Conversely, premature removal of drains can cause fluid buildup, which hinders pulmonary function. **Continuous pH and impedance monitoring of chest drainage can facilitate quantitatively informed and prudent decision-making when determining the optimal time for tube removal.**



**Figure 1.** Schematic showing inline attachment of Origin™ to traditional chest drainage systems.

## Integration of Origin™ with Traditional Chest Tubes




Impedance sensors can be used to objectively monitor and track postoperative air leaks<sup>4</sup>. reducing the risk of prolonged air leaks pH sensors can be used to infer the inflammatory processes associated with empyema and other postoperative cardiopulmonary complications<sup>5</sup>.


FluidAI Medical developed a non-invasive, sensor-based platform (Origin™) that can be modified to connect in-line to traditional/analog chest drainage systems for real-time 24-hour monitoring of drained effluent. **This study assesses the ability of Origin™ to accurately measure the pH of pleural fluid compared with the widely used Radiometer ABL800 Flex blood gas analyzer (BGA).**

## Methods


Thoracotomy was performed on 9 Yorkshire pigs (Mean Age: 4 Months, Mean Weight : 31 kg)



Pleural fluid was collected from the left and right pleural spaces.



The pH of each sample was analyzed using Origin™ and BGA.



Pearson's correlation, linear regression analysis, and Bland-Altman (BA) analysis were used to compute inter-technique agreement.



## Results

The Pearson coefficient ( $r = 0.965$ ,  $p < 0.05$ ) and the coefficient of determination ( $r^2 = 0.931$ ,  $p < 0.05$ ) indicated a strong linear correlation between pH data collected using the two systems within a measuring range of 6.000-8.200. The linear regression model had a slope of 0.820 and an intercept of 1.037.

BA analysis revealed that the mean estimated bias  $\pm$  standard deviation between pH measurements obtained using Origin™ and BGA was  $0.382 \pm 0.084$ , (95% CI: 0.293 – 0.470). The reported upper and lower limits of agreement were 0.547 (95% CI: 0.386 – 0.708) and 0.216 (95% CI: 0.055 – 0.377) respectively.

## Conclusion

**Origin™ can seamlessly integrate with existing traditional/analog chest drainage systems. pH measurements made using Origin™ strongly correlate with BGA measurements.**

Continuous pH and impedance analysis of pleural effluent using Origin™ can help improve the management of traditional chest tube drainage systems.

## References

- [1] Utter GH. The rate of pleural fluid drainage as a criterion for the timing of chest tube removal: theoretical and practical considerations. *Ann Thorac Surg.* 2013;96(6):2262-2267. doi:10.1016/j.athoracsur.2013.07.055
- [2] Anderson D, Chen SA, Godoy LA, Brown LM, Cooke DT. Comprehensive Review of Chest Tube Management: A Review. *JAMA Surg.* 2022;157(3):269-274. doi:10.1001/jamasurg.2021.7050
- [3] Kwiat M, Tarbox A, Seamon MJ, et al. Thoracostomy tubes: A comprehensive review of complications and related topics. *Int J Crit Illn Inj Sci.* 2014;4(2):143-155. doi:10.4103/2229-5151.134182
- [4] DeArmond DT, Das NA, Restrepo CS, et al. Intrapleural Impedance Sensor Real-Time Tracking of Pneumothorax in a Porcine Model of Air Leak. *Semin Thorac Cardiovasc Surg.* 2020;32(2):357-366. doi:10.1053/j.semtcvs.2019.10.004
- [5] Houston MC. Pleural fluid pH: diagnostic, therapeutic, and prognostic value. *Am J Surg.* 1987;154(3):333-337. doi:10.1016/0002-9610(89)90623-5